

PLATINUM



Description

Single phase Hydrogen Peroxide solution for powerful cleaning, disinfection, neutralisation, protein removal and storage of contact lenses.

Pack size available

60ml, 250ml (FB), 360ml

Key features

- The formulation deeply cleans, loosens debris and deposits
- Suitable for use with all soft and RGP contact lenses, less suitable for coloured lenses
- Good for those with allergic responses and sensitive eye
- Platinum storage case ensures initial time for deep cleaning in peroxide before the neutralising process.

Ingredients

Hydrogen Peroxide 3%, Dibasic sodium phosphate, Monobasic sodium phosphate, stabilized with Ethylene diamine tetra acetic acid 0.03% and Phosphonic acid, Natrium chloride, Boric acid, Purified water.



A highly effective disinfecting solution to help keep your contact lenses clean



Formulated specially to prevent deposits



Preservative free and ideal for those with more sensitive eyes



Good for allergies



Suitable for all contact lenses, soft and RGP

Comparing your old Soleko XXX to your new Soleko Platinum

		Soleko XXX	Soleko Platinum
Unique selling point		• Single phase hydrogen peroxide solution (no tablet) • Powerful disinfection cleaning system for all types of lenses • Good for sensitive eyes	• Single phase hydrogen peroxide solution (no tablet) • Powerful disinfection cleaning system for all types of lenses • Good for sensitive eyes
Formulation	Disinfecting ingredient	Hydrogen peroxide 3%	Hydrogen peroxide 3%
	Cleaning ingredients	Poloxamer	-
	Neutralisation method	Disc in case	Disc in case
	Other	NaCl, water, phosphates HEDP (stabiliser and chelating)	Phosphate buffer, EDTA, NaCl, boric acid, water
Soaking time		At least 6 hrs or overnight	At least 6 hrs or overnight
Storing period		7 days	24 hrs
Discard time after open		90 days	60 days